ESOTEST® MULTI

SAFETY AND PERFORMANCE INFORMATION

(EN)

1 INDICATIONS FOR USE

The Esotest® Multi Temperature Monitoring System is composed of the Esotest® Multi temperature monitor and the Esotherm® Multi temperature probe and is intended for the continuous detection, measurement and visualization (in °C) of esophageal temperature. The intended environments of use are operating rooms and interventional electrophysiology rooms.

1.1 Essential performances

The combination of Esotest® Multi and Esotherm® Multi, properly interconnected, constitutes a clinical thermometer for continuous use, therefore the essential performance of the system consists in measuring and displaying the temperatures of the probe sensors with a predefined accuracy and in alerting the user if for any reason such measurement is impossible. In particular:

- In absence of hardware faults or significant disturbances, the system ensures the continuous measurement and visualization of the temperatures detected by the sensors, with the accuracy declared in the performance specification summary (see section 4).
- If the system is not able to perform the measurement with the expected accuracy (for example due to hardware faults or high intensity environmental disturbances), the temperature visualization stops and the event is notified to the user with appropriate warning signals.

With the optional accessory Esotest® Multi adapter cable, Esotest® Multi can be connected to an external mapping system for visualizing the position of the probe sensors inside the patient's body.

WARNING!

The technology used by Esotest® Multi for the measuring process requires that the handle of the probe and that of the patient cable are approximately at the same temperature. Hence, it is recommended to wait at least 15 minutes between the instant in which the components are brought in the working environment and the start of a measurement session, especially if the probe or the patient cable have long stationed in environments with temperature and/or humidity significantly different from those of use.

For the same reason, it is recommended not to place the two handles in such a way that during the measurements only one of the two objects remains for a prolonged time in proximity to or in contact with a heat or cold source. Failure to observe these precautions can lead to an accuracy of the instrument lower than the one stated in the technical specifications.

WARNING!

Monitoring the esophageal temperature does not guarantee the safety of the heating/cooling process under observation. It is therefore essential that the devices responsible for the heating/cooling process are always operated according to the safety recommendations listed in their instruction manuals.

1.2 Contraindications and side effects

There are no specific contraindications or side effects related to the use of Esotest® Multi. However, the insertion of the probe by an unskilled operator can cause discomfort and slight abrasions to the patient. To avoid that risk, it is essential that Esotest® Multi is used only by expertised medical personnel.

The insertion of the Esotherm® Multi probe in the esophagus is not recommended if the esophagus is deformed or irritated.

2 DEVICE INSTALLATION

2.1 Package contents

The **Esotest® Multi** package contains:

- 1 Main unit
- 1 Patient cable
- 1 Equipotential cable
- 1 External power supply
- 1 AC power cord
- 1 Printed copy of the Instructions for Use
- 1 Test certificate
- 1 Warranty certificate

When opening the package for the first time, please check that:

- The wrapping is intact
- All accompanying documents are present
- The device has been regularly tested and it is fully compliant (as certified by the test certificate)
- Cables and accessories are all present and undamaged
- The AC power cord is suitable for the country of use

The **Esotherm® Multi** package contains a sterile double pouch with the following components inside:

- 1 Esotherm® Multi probe equipped with one Ø 0,35 mm stylet (dark blue handle) already inserted
- 1 spare Ø 0,40 mm stylet (yellow handle) placed inside a plastic ring
- 1 Printed copy of the Instructions for Use

The **Esotest® Multi Adapter Cable** package contains:

- 1 Esotest® Multi Adapter Cable

In case of damage/tampering of the packaging or in the absence/damage of one of the components listed above, it is recommended to contact the customer service immediately.

2.2 Assembling the Esotest® Multi monitoring system

Figure 1 shows all the components of Esotest® Multi and how they must be interconnected before switching on the device.


Figure 1: Components of Esotest® Multi

Warnings pertaining to the whole system:

**WARNING!**

When interconnecting the various parts of the machine it is necessary to verify that none of the cables (AC and DC power cables, equipotential cable, patient cable, adapter cable) appears visibly damaged. The presence of exposed copper wires or tears in the cable sheathing can compromise the safety and the proper operation of the system.

**WARNING!**

The medical grade external power supply is an integral part of Esotest® Multi and it's critical to ensure the safety of the system. Consequently, it is recommended to power Esotest® Multi exclusively with the power supply included in the package.

**WARNING!**

To avoid the risk of electric shock, Esotest® Multi must be connected only to supply mains with protective earth. Only in this way the external power supply can ensure the safety level required by medical standards.

**WARNING!**

Some of the insulation barriers that guarantee the safety of the patient cable are implemented inside the main unit. Consequently, it is recommended to keep the plug of the patient cable inserted in its socket (A5) as long as the temperature sensors touch the patient's body.

**WARNING!**

The circuits inside the handle of the patient cable are essential to ensure the safety and the accuracy of Esotest® Multi. Consequently, it is recommended not to subject the handle to impact or other mechanical shocks, nor to tamper with it in any way.

**WARNING!**

During the operation of Esotest® Multi the sensors (B2) of the temperature probe touch inner parts of the patient which must never be exposed to uncontrolled electrical signals. For this reason, before inserting the probe sensors inside the human body it's necessary to verify the integrity of the patient cable (C), taking special care to protect socket C2 with its rubber stopper if that connector is not in use. For the same reason, never disconnect the patient cable from the main unit until the temperature sensors are extracted from the patient or at least the temperature probe has been detached from socket C1.

Warnings pertaining to the temperature probe:

**WARNING!**

The sterility of Esoterm® Multi may be compromised if it is stored incorrectly. Before opening the package, check that it has not been damaged and that its sterility has not been compromised. Do not use products contained in broken or damaged packages. Do not use sterile products after their expiration date (as indicated on the packaging).

**WARNING!**

Check the correct insertion of the probe into the patient cable connector and the effectiveness of the connector retention by pulling the connector of the probe itself. The improper handling (traction, torsion, bending, etc.) of the lead body and in particular of the distal segment, can cause damage and loss of functionality.

**WARNING!**

If the probe is lubricated to facilitate the insertion of the sensors into the esophagus, it is essential to prevent the lubricant from coming into contact with the metal terminals of the probe plug (B1). The penetration of liquids inside the socket reserved for the probe may compromise the operation of the measuring system and damage the electronic circuits of the patient cable.

**WARNING!**

When the temperature sensors are inside the patient's esophagus, the stylet, the plug (B1) on the handle and the nearby metal hooks provide an electrical path towards the internal organs of the human body. To ensure patient safety, it is therefore recommended to take all the necessary precautions (for example, preventively touching the equipotential node of the operating room) to avoid an electrostatic discharge between the hands of the operator and the metal parts of the probe.

**WARNING!**

The use of components and/or accessories other than the supplied parts during probe insertion may damage the probe and/or degrade its performance

**WARNING!**

Do not submit to MRI a Patient with Esotherm® Multi inserted.

Warnings pertaining to the adapter cable:

**WARNING!**

The contacts inside the adapter cable socket provide unprotected electrical paths toward the probe electrodes. For this reason, after the introduction of the electrodes into the human body it is essential to avoid sending through the adapter cable socket electrical signals hazardous for the patient's health, such as electrostatic discharges or signals produced by non-medical devices. For the same reason, after each use it is recommended to cover the socket with its rubber plug.

**WARNING!**

When Esotest® Multi is connected to an external system, the safety and the performance of the device are guaranteed only if the electromagnetic emissions of that system comply with standard IEC 60601-1-2 and its ports are classified as a CF-type, defibrillator-proof applied parts.

**WARNING!**

An external device connected to the socket for mapping systems via the adapter cable may interfere with the measurement circuits inside the patient cable. If during the connection one or both devices exhibit abnormal behavior, it is recommended to interrupt the temperature measurement as long as the external device is operating.

**WARNING!**

When the adapter cable is inserted in the handle of the patient cable, the single pole connectors provide unprotected electrical pathways toward the probe electrodes, so it's required to plug them into the input ports of the mapping system before inserting the probe in the patient's esophagus and, in general, to keep the plugs away from any source of electrical signals that may be harmful to the patient.

3 PLACING ESOTEST® MULTI IN THE WORKING ENVIRONMENT

The following sections describe the preventive measures required to ensure that the device, when placed in its working environment, always operates safely and at the level of performance stated in its technical specifications.

3.1 Emissions and immunity to electromagnetic interference

Esotest® Multi does not emit any kind of ionizing radiation. With regard to non-ionizing radiations, the device complies with the current regulations relating to electromagnetic compatibility (EMC). The emissions of the device are in line with the standards for non-domestic working environments and should not compromise the

performance of surrounding electronic instruments.

If Esotest® Multi is used in an environment characterized by electromagnetic disturbances comparable to those defined by the IEC 60601-1-2 standard, the device is able to guarantee the essential performances described in section 1.1. In presence of disturbances of greater magnitude or with characteristics that the standard does not take into consideration (for instance, disturbances that penetrate by conduction through the electrodes of the probe), Esotest® Multi may behave abnormally. For example:

- The system may detect temperature changes characterized by a scale and rapidity incompatible with the processes under observation
- It may occur a contemporary and sudden increase or decrease of all the measured temperatures for no apparent reason
- With all the elements of the system properly interconnected, a dialog window may suddenly appear, displaying a message which notifies the interruption of the communication with the patient cable
- The touch-screen may stop working

In particular, Esotest® Multi may be affected by electrostatic discharge (ESD). In an environment likely to cause ESD, such as a room with carpeted floor, the operator must take care to get rid of static electricity on his/her body before touching the device. It is recommended that attending health care professionals discharge any static electricity by touching a large metal or conductive, grounded surface before touching the patient, the cable, the leads, or the Esotest® Multi.

Additional disturbances not taken into consideration by the standard may be generated/focused inside the patient's body and mainly involve the probe sensors, or they may be diffused in the environment by equipments placed near the main unit of the device.

- With regard to the sources of interference inside the patient's body, the device is immune to electromagnetic disturbances generated by radiofrequency ablaters, navigation systems and cardiac stimulators. If the patient is subject to a defibrillation discharge, the measured temperatures may show an extremely short positive or negative peak at the instant of the discharge. There is no data concerning the behavior of Esotest® Multi in presence of additional sources of electromagnetic disturbances located/focused in the patient's body, such as electrocauteries, RFIDs or diathermy machines, so it's not recommended to perform temperature measurements while these devices are active
- With regard to sources of disturbance outside the operating table, the system may behave abnormally in close proximity to electric motors, power lines, or other electronic equipment in operation, especially portable radio communication devices. In such cases, it is recommended to identify the source of the interference and to turn it off or move it away.

As a rule, once the source of the disturbance is removed, Esotest® Multi resumes operation without requiring further user intervention. If, on the other hand, the disturbance has caused the deactivation of the touch-screen or the interruption of the communication with the patient cable, the user can manually restart both functions as follows:

- For restoring the touch-screen functionality, press the on/off button twice in succession
- For restoring the communication with the patient cable, first check whether the patient cable is correctly connected to the main unit of Esotest® Multi. If the connection does not have any issue, follow the instructions in the message window.

If the abnormal behaviors described above persist even after the removal of all possible causes of interference, contact the customer service.

3.2 Use in a system composed of multiple devices

Esotest® Multi is commonly used in the context of a system composed of multiple devices. In this case, assessing the conformity of the system to specific electrical safety and electromagnetic compatibility requirements is responsibility of the system designer. As a rule, it is recommended to connect all devices in the same workplace to the equipotential node of the electrical system. If Esotest® Multi is employed together with other devices with applied parts, it is essential that these applied parts shall be of CF type, so as to avoid a potentially dangerous increase of leakage currents on the patient. More generally, the use of

equipments not complying with safety requirements equivalent to those of Esotest® Multi can lead to a reduced level of safety of the resulting system. Additionally, in order to minimize the possible degradation in the performance of Esotest® Multi, it is recommended to keep all power and signal cables of the device away from other power and signal cables which may be present in the environment.

4 PERFORMANCE SPECIFICATION SUMMARY

Main unit and reusable components:

Device definition	<p><i>Model:</i> Esotest® Multi <i>Description:</i> clinical programmable temperature monitor <i>Class and type (in accordance with EN 60601-1 standard):</i> Class I mains powered equipment, transportable <i>Type of applied part:</i> CF, defibrillation-proof <i>Conditions of use:</i> device for continuous use <i>Usage environment:</i> operating rooms and interventional electrophysiology rooms</p>
Temperatures measurement	<p><i>Number of temperature sensors:</i> 3, 5 or 7 <i>Temperature measurement range:</i> from -15°C to 75°C <i>Measurement accuracy:</i> ±0.3°C from 10°C to 45 °C (rated range) ±0.4°C from 0°C to 75°C (extended range) ±0.4°C from -15°C and 0°C is, verified by the Manufacturer <i>Resolution:</i> 0.1°C <i>Response Time*:</i> 1s <i>Heating/Cooling Transient Time (according to ISO_80601-2-56):</i> 1.8s</p> <p>* The "Response Time" or "Time Constant" is defined as the time taken by the system to detect 63.2% of an instantaneous temperature change</p>
Temperature recording	<p><i>Maximum duration of a single recording:</i> 8 hours <i>Sampling rate in archived treatments:</i> 1 sample/second for each sensor</p>
Main alarm thresholds	<p><i>Range of the lower alarm threshold (Tmin):</i> from 12°C to 24°C in steps of 1°C <i>Range of the upper alarm threshold (Tmax):</i> from 37°C to 41°C in steps of 0.5°C <i>Tmin default value:</i> 18°C <i>Tmax default value:</i> 39°C</p>
Operating conditions	<p><i>Room temperature:</i> from +10°C to +40°C <i>Relative humidity, uncondensed:</i> from 30% to 75% <i>Atmospheric pressure:</i> from 700 mbar to 1060 mbar</p>
Storage and transport conditions	<p><i>Room temperature:</i> from 0°C to +50°C <i>Relative humidity, uncondensed:</i> up to 80% <i>Atmospheric pressure:</i> from 700 mbar to 1060 mbar</p>

Temperature probe:

Models, dimensions and number of sensors	<p><i>Model:</i> Esotherm® Multi <i>Total esophageal tube length:</i> 83 cm <i>Usable length:</i> 75 cm <i>Esophageal tube diameter:</i> 2.3 mm / 7 French <i>Maximum electrode diameter:</i> 3.6 mm / 11 French <i>Number of temperature sensors:</i> 3 (Ref. 2613535H) 5 (Ref. 2613558H) 7 (Ref. 2613576H) <i>Distance between electrodes / Measurement array length:</i> 5 mm / 25 mm (3 sensors) 8 mm / 47 mm (5 sensors) 8 mm / 69 mm (7 sensors)</p>
Storage and transport conditions	<p><i>Room temperature:</i> from 0°C to +50°C <i>Relative humidity, uncondensed:</i> from 10% to 80%</p>
Expiry date	48 months after the sterilization date

Rev.	Revision Status / Description	Date
A	First Issue	28/09/2022
B	The document layout has been changed to conform to model 99501185	02/08/2023

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