REUSABLE NEUTRAL ELECTRODES FOR ELECTROSURGERY INFORMATION FOR SAFETY AND PERFORMANCE (EN)

1. PRODUCT DESCRIPTION

Reusable neutral electrodes for electrosurgery are non-sterile medical devices, class IIb according to EU Regulation 2017/745 (MDR), intended to provide an electrical return path for the high-frequency current introduced in the patient by the generator through the active part in electrosurgical procedures.

2. INDICATIONS AND USERS

Reusable neutral electrodes are medical devices with a relatively large surface area applied to the patient's body, so as to reduce the density of the current returning from the generator (ESU - ElectroSurgical Unit) during electrosurgical procedures and thus avoid undesirable effects such as burns. FIAB's electrodes can be used also on multiple patients, upon adequate cleaning according to the instructions for use.

Neutral electrodes must be applied to the patient, connected to the electrosurgical unit and used only by medical personnel qualified for such procedures, i.e., physicians specialised in electrosurgery or operating room nurses.

3. INTENDED PURPOSE

Used in connection with high-frequency generator during electrosurgical procedures, the neutral electrode provides an electrical return path for the high-frequency current introduced in the patient through the active part.

4. ESSENTIAL PERFORMANCE

Provide the preferential return path for the circuit composed of the electrosurgical unit, pencil and electrode, preventing dangerous dissipation of electrical energy to other conductive parts.

5. CONTRAINDICATIONS

There are no particular contraindications to the use of neutral reusable electrodes.

6. POSSIBLE COMPLICATIONS

Possible causes of burns include, but are not limited to:

- Excessive current density in the patient's tissue,
- Incorrect placement of the electrode,
- Accidental contact of the patient with electrically conductive parts,
- Direct contact between the cables and the patient's skin with resulting capacitive effect,
- Combustion of inflammable disinfectant agents,
- Combustion of inflammable narcotic gases.

6.1 Residual risk

After risk control measures have been taken, there are no residual risks to the use of the device.

7. PRECAUTIONS AND WARNINGS

To reduce the complications listed above the following precautions need to be taken.

7.1 Packaging

WARNING!

Check the integrity of the electrode packaging: do not use the product in case of damages or visible defects.

WARNING!

The product must be stored in its original packaging under the environmental conditions (temperature and relative humidity) specified on the label.

WARNING!

Do **NOT** use in case of damaged packaging.

7.2 Equipment and accessories

WARNING!

Connect the electrode to the generator using the connector cable.

WARNING!

For connections and correct functioning, follow the instructions for use supplied by the manufacturer of the ESU unit and pencil being used.

WARNING!

High-frequency electrosurgical generators, pencils and accessories to be used with this product must comply with the regulations in force.

WARNING!

Check that the generator alarm system is working before starting the procedure.

WARNING!

Pacemakers or other active implants can be damaged by an electrosurgical current.

7.3 Electrical and electromagnetic safety

WARNING!

The current must never flow in a transversal direction across the body, nor cross the chest.

WARNING!

Do not submit a patient with implanted pacemaker or other active implant to an electrosurgical current without first consulting a cardiologist or another qualified person.

EMC DISTURBANCES – Refer to the generator operating instructions for information about:

- radio interferences which may be produced during operation and disrupt nearby equipment,
- information on sensitivity to external electromagnetic disturbances that could cause degradation of the desired surgical effects.

7.4 Performance and use

WARNING!

Choose a well-vascularised portion of skin close to the area to be operated: do not place the electrode over scars or close to metallic prostheses or near ECG electrodes.

WARNING!

Shave, degrease and dry the chosen surface.

WARNING!

Make sure that the entire surface of the electrode is in contact with the patient's skin.

WARNING!

During the procedure, always select the lowest energy level possible. If the coagulation ability of the pencil electrode is lower than usual, do not increase the high-frequency output before checking:

- the correct positioning of the neutral electrode (perfect adherence on the entire surface),
- the correct connection of cables and their connectors,
- the correct activation of the starting keys (hand-switch or foot-control),
- that there is no damage to the insulation of the cables,
- that the pencil electrode is not dirty.

WARNING!

Use is not recommended if high-power values are applied for long times and in presence of electrolytic solutions.

WARNING!

If the patient is "re-positioned", check the contact of the electrode with the patient's skin and the cable connections.

WARNING!

Incorrect positioning of the electrode can cause burns: follow these instructions.

WARNING!

Do not apply the electrode over wounds or scars, close to metallic prostheses or in the vicinity of ECG electrodes and in areas of potential liquid flow. The patient must be positioned on a dry and electrically insulated surface. The patient must be kept insulated from conductive parts, and the operation table must be suitably "grounded". Use dry gauze to avoid areas of skin coming in contact with each other.

WARNING!

FIAB recommends checking the integrity of the device before each use and to avoid exceeding the limit reported above even if there is no visible damage.

WARNING!

Following the cleaning procedure, FIAB guarantees exclusively the device functionality after cleaning. The efficacy of cleaning itself is dependent on the final user.

WARNING!

DO NOT STERILIZE.

WARNING!

The product may be contaminated after use, constituting a biological hazard. The device and its parts must be handled according to medical procedure in order to prevent injury or infection.

Rev.	Revision Status / Description	Date
A	First Issue	31/07/2023

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