# **MASKS FOR OXYGEN AND AEROSOL THERAPY**

## **INFORMATION FOR SAFETY AND PERFORMANCE**

### **(EN)**

## **INDICATIONS**

Therapeutic administration of oxygen.  
Disposable product (for a single patient), to be used continuously for a period not exceeding 30 days.  
The oxygen source must be in according to the regulations in force.  
The products must be used by qualified healthcare personnel.

## **CONTRA-INDICATIONS**

No contraindications for oxygen therapy are described.

## **PRECAUTIONS AND WARNINGS**

Reuse of the product may result in:

- 1 presence of biological residues that might cause cross-infections,
- 2 alteration of materials,
- 3 loss of initial functional features of product.

Before connecting to the oxygen source, check its operation.

Accidental crushing or kneeling of the connection tubing do not allow oxygen to flow properly and may cause the disjunction of the connector. It is therefore recommended that the users control carefully.  
The same attention must be paid even to the models which are supplied with a sure-flow tubing for which the risk of disjunction exists even though in a reduced way.

In case of product containing phthalates (DEHP). Even if risks to human health are not confirmed, the use of the devices on pregnant/nursing or child should be evaluated by your doctor.

The concentrated oxygen can ignite combustible materials. To avoid risk of fire or explosion keep away from possible sources of ignition, and refer to the precautions for use of the equipment used for the delivery of oxygen.

Rev.	Revision Status / Description	Date
A	First Issued	08/05/2023
B	Introduction of form 99501185	03/07/2023
C	Updated with the header corrected	17/10/2023

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